

# High Frequency Power with the ease of "touch" Anatomical programming!

The power of the InnoVet Select™ High Frequency X-ray unit now includes a new option. Introducing the new anatomically-programmed InnoVet Select High Frequency. The High Frequency generator that is conveniently mounted above the Select's SafeRay™ collimator.



This new anatomically-controlled generator provides the InnoVet operator maximum ease-of-use to make imaging with the InnoVet Select HF fast and simple.

## It's as easy as 1-2-3

- Step 1:** From the console, select the anatomical region of your choice.
- Step 2:** Then, select the particular view.
- Step 3:** Select the body part thickness.



Now you're ready to produce a good radiograph. The software for the anatomically-controlled generator was designed for InnoVet by board-certified Veterinary Radiologists. The pre-programmed techniques have been proven through hours of field testing. In addition, this new system allows for the adjustment of techniques for different film speeds or density preferences.

### Anatomical Generator Specifications:

- Frequency rating: 20 kHz
- Power Rating: 30kW (300mA @ 100 kVp)
- kVp Range: 40-125 in 1 kV steps
- mA Range: 50 - 300 (50S, 100S, 100L, 300L)
- mAs Range: 0.2 - 400
- Adjustable exposure signal duration
- Number of procedure categories: 8
  - Featuring:** thorax, abdomen, two extremity settings, head/spine, shoulder/pelvis, exotics, and special procedures
- 50 commonly used, pre-programmed positions to select from
- Over 1,000 pre-programmed techniques



## InnoVet™

Look to the Paw™

[www.innovet4vets.com](http://www.innovet4vets.com)

1.800.729.9729